



LarySeal™

Comprehensive Range of Laryngeal Mask Airways



FEDIMED
Sociedad Anónima

LARYSEAL MULTIPLE

LARYSEAL BLUE

LARYSEAL CLEAR

LARYSEAL MRI

LarySeal™ Multiple

A silicone Reusable Laryngeal Mask Airway designed to be sterilized and used up to a recommended 40 times. The LarySeal Multiple's cuff adapts to the contour of the oropharyngeal area to provide a secure seal.

- Medical Grade Silicone
- Autoclavable

Colour Coded Pilot Balloon

Ease of identification of sizing by colour coded pilot balloons.

Secure Seal for Airway Management

Highly Elastic Silicone material allows the cuff of the mask to form itself to the contour of the oropharyngeal area and provide a seal.

Universal Connection


15mm M connector provides a universal connection.

Crush Resistant

The main tube wall thickness reduces the risk of crushing.

Easy indication of the cuff pressure through pilot balloon

Thin wall pilot balloon indicates the inflation of the cuff.



Autoclavable, Reusable

Durable high temperature resistant material. Recommended for up to 40 times autoclaving.

PRODUCT CODE	DESCRIPTION	PACK SIZE
038-94-110	LarySeal Multiple size 1	10
038-94-115	LarySeal Multiple size 1.5	10
038-94-120	LarySeal Multiple size 2	10
038-94-125	LarySeal Multiple size 2.5	10
038-94-130	LarySeal Multiple size 3	10
038-94-140	LarySeal Multiple size 4	10
038-94-150	LarySeal Multiple size 5	10

LarySeal™ Blue

The first Silicone Single Patient Use Laryngeal Mask Airway introduced to provide the benefits of the well accepted reusable LarySeal with the added benefits of single use, protecting against cross contamination.

- Medical Grade Silicone
- Single Use
- Cost Effective

Colour Coded Pilot Balloon

Ease of identification of sizing by colour coded pilot balloons.

Less risk of slipping or dislocation

The silicone matt surface finish of the cuff provides a good adhesion to the oropharyngeal area and reduces the risk of movement and dislocation of the cuff.

Universal Connection


15mm M connector provides a universal connection.

Crush Resistant

The main tube wall thickness reduces the risk of crushing.

Easy indication of the cuff pressure through pilot balloon

Thin wall pilot balloon indicates the inflation of the cuff.



Reduce the risk of cross contamination

Sterile product in secure packaging.

PRODUCT CODE	DESCRIPTION	PACK SIZE
038-94-210	LarySeal Blue size 1	10
038-94-215	LarySeal Blue size 1.5	10
038-94-220	LarySeal Blue size 2	10
038-94-225	LarySeal Blue size 2.5	10
038-94-230	LarySeal Blue size 3	10
038-94-240	LarySeal Blue size 4	10
038-94-250	LarySeal Blue size 5	10

LarySeal™ Clear

LarySeal Clear is a single patient use PVC Laryngeal Mask Airway, engineered using extra soft matt textured materials to replicate as closely as possible that of silicone, providing Anaesthetists the same feel as LarySeal Blue and LarySeal Multiple.

- Medical Grade PVC
- Single Use
- Cost Effective

Crush Resistant

The main tube wall thickness reduces the risk of crushing.

Colour Coded Pilot Balloon

Ease of identification of sizing by colour coded pilot balloons.

Less risk of slipping or dislocation

The satin surface finish of the cuff provides a similar texture to silicone, thus providing ease of insertion while providing a secure seal.

Universal Connection


15mm M connector provides a universal connection.

Easy indication of the cuff pressure through pilot balloon

Thin wall pilot balloon indicates the inflation of the cuff.

Cost Effective

Eliminate the need for the costly and time consuming process of cleaning and re-sterilizing.



LarySeal MRI

Providing the same benefits of LarySeal Clear with the added benefit of a complete metal free inflation valve. LarySeal MRI can be used within any MRI suite and will not distort the images.

PRODUCT CODE	DESCRIPTION	PACK SIZE	PRODUCT CODE	DESCRIPTION	PACK SIZE
038-94-310	LarySeal Clear size 1	10	038-94-410	LarySeal MRI size 1	10
038-94-315	LarySeal Clear size 1.5	10	038-94-415	LarySeal MRI size 1.5	10
038-94-320	LarySeal Clear size 2	10	038-94-420	LarySeal MRI size 2	10
038-94-325	LarySeal Clear size 2.5	10	038-94-425	LarySeal MRI size 2.5	10
038-94-330	LarySeal Clear size 3	10	038-94-430	LarySeal MRI size 3	10
038-94-340	LarySeal Clear size 4	10	038-94-440	LarySeal MRI size 4	10
038-94-350	LarySeal Clear size 5	10	038-94-450	LarySeal MRI size 5	10

LarySeal Features


The LarySeal range of Laryngeal Mask Airways has the following additional common benefits:

Reduced risk of blockage




The deep design of the patient port reduces the risk of blockages by soft tissues. The open patient end port allows the use of endoscope without restriction.

Full glottis view



The supraglottic airway device allows for a clear view of the glottis using an endoscope.


Printed info on the main tube



Guidance for appropriate air volume and size is printed on the main tube, allowing for quick and easy reference.

Clearly labelled packaging

Instructions for use are clearly shown on the packaging in both text and diagrammatic form.



LarySeal in position

The diagram shows the LarySeal in position and demonstrates how the soft material of the cuff forms itself to the oropharyngeal area.

Colour coded LarySeal Cuff Inflator

The LarySeal Cuff Inflator has clearly marked increments, each one specifically colour coded to match the size colour coding of the LarySeal pilot balloon. Thus enabling easier and quicker identification of the correct measure for inflating the cuff of the LarySeal.




Laryngeal Mask Airway Colour Coding

All of the LarySeal range has a colour coded pilot balloon which corresponds to the size of the Laryngeal Mask Airways, providing a clear visual indication of the size, therefore facilitating quicker selection.

Colour Coded Pilot Balloon

Ease of identification of sizing by colour coded pilot balloons.



THE LARYSEAL PILOT BALLOON COLOUR CODING SYSTEM

	White Pilot Balloon - 1	 <4ml/60cm H ₂ O	 Size 1 < 5Kg
	Blue Pilot Balloon - 1.5	 <7ml/60cm H ₂ O	 Size 1.5 5Kg - 10Kg
	Light Green Pilot Balloon - 2	 <10ml/60cm H ₂ O	 Size 2 10Kg - 20Kg
	Green Pilot Balloon - 2.5	 <14ml/60cm H ₂ O	 Size 2.5 20Kg - 30Kg
	Orange Pilot Balloon - 3	 <20ml/60cm H ₂ O	 Size 3 30Kg - 50Kg
	Red Pilot Balloon - 4	 <30ml/60cm H ₂ O	 Size 4 50Kg - 70Kg
	Yellow Pilot Balloon - 5	 <40ml/60cm H ₂ O	 Size 5 70Kg - 100Kg

The LarySeal Cuff Inflator has been specifically designed for use only with the LarySeal range of Laryngeal Mask Airways, indicating the maximum amount of air by individual colour coded sizes.

	PRODUCT CODE	DESCRIPTION	PACK SIZE
	038-94-001	LarySealColor Coded Cuff Inflator	10

Details and specifications contained in this brochure are correct at time of printing (May 2007), but this publication should not be regarded as an infallible guide.

For further information please contact:

Fedimed S.A.

Av. Ricardo Balbín 3402 (C1430AAS) - Buenos Aires - Argentina

Tel. 4543-4114 | FAX 4544-0008

www.fedimed.com.ar | info@fedimed.com.ar