

LUNAR - the low cost reusable laryngoscope system.



Bright, focused replaceable bulb

A powerful replaceable bulb provides excellent illumination which is directed towards the airway.

Stainless Steel Blade

Produced to identical profiles as our Orion reusable system. Fully autoclaveable.

Standard fitting plastic block

The high impact block allows low cost production without compromising rigidity or strength. It is also marked with the blade profile and size. Made to popular ISO 7376/1 standard light fittings, interchangeable with all quality standard light systems.

3 Ball Bearings

Provide a secure fitting with a reassuring click

Lunar - the new reusable standard light laryngoscope system from Timesco providing excellent quality and superb clinical performance at an unbeatably low price.

As one of the world's leading laryngoscope brands we recognised the need for an economically priced system for emergency and pre-hospital applications.

The Lunar system has been designed to be indistinguishable in use to our full metal systems but the economy of the plastic block allows a much lower production cost and selling price. Less than half the cost of a comparable reusable system!

The plastic block is produced from a special

polymer which provides maximum strength which will not deform or snap under pressure. All Lunar blades and handles are made to the popular ISO 7376/1 standard light fitting and are fully autoclaveable to 50 cycles at 134°C for 3 minutes. Sterilisation may also be performed by Ethelne Oxide or cleaning with alcohol.



2904.150.10 Lunar Macintosh Laryngoscope Blade 1 (92mm x 14mm)
2904.150.15 Lunar Macintosh Laryngoscope Blade 2 (100mm x 17mm)
2904.150.20 Lunar Macintosh Laryngoscope Blade 3 (130mm x 18mm)
2904.150.25 Lunar Macintosh Laryngoscope Blade 4 (154mm x 18mm)

2904.185.05 Lunar Miller Laryngoscope Blade 0 (78mm x 9.5mm)
2904.185.10 Lunar Miller Laryngoscope Blade 1 (105mm x 9.5mm)
2904.185.15 Lunar Miller Laryngoscope Blade 2 (158mm x 13mm)
2904.185.20 Lunar Miller Laryngoscope Blade 3 (198mm x 13mm)
2904.185.25 Lunar Miller Laryngoscope Blade 4 (207mm x 14mm)

2904.300.05 Lunar Laryngoscope Handle Size Small AA
2904.300.10 Lunar Laryngoscope Handle Size Medium C

2905.325.10 Laryngoscope Bulbs Small 2.5V Tungsten Packet of 6
2905.325.20 Laryngoscope Bulbs 2.5V Standard Tungsten Packet of 6
(Small bulbs are for use with Miller sizes 0 & 1 blade only).

2904.900.50 Lunar Laryngoscope Set 1 Mac 1, 2, 3 Blades, Medium Handle & Hard Plastic Case
2904.900.75 Lunar Laryngoscope Set 2 Mac 2, 3, 4 Blades, Medium Handle & Hard Plastic Case



Timesco
SURGICAL + MEDICAL
INSTRUMENTS FOR LIFE



SURGICAL • LARYNGOSCOPY • FOOTCARE • EMERGENCY • DIAGNOSTIC • MIDWIFERY • DENTAL

Timesco of London Ltd, Timesco House, 1 Knights Road, London E16 2AT
Telephone: +44 (020) 7511 9960, Fax: +44 (020) 7511 7888, E-Mail info@timesco.com
www.timesco.com

Copyright April 2004, Timesco, Instruments for Life, Lunar are registered trademarks of Timesco of London Ltd.



View Timesco's
range of products
on the World Wide
Web.
www.timesco.com