



THE NEW OPTIMA XL XENON HANDLE OVER 100% MORE LIGHT!

Our new Optima XL laryngoscope handle provides the same excellent clinical performance as our Optima handle, but with a brighter, whiter illumination.

- *Over 100% increase in light*
- *Whiter illumination*
- *Powerful Xenon bulb*
- *Simple clip in bulb fitting*
- *Fully compatible with Timesco Optima laryngoscope blades and other ISO 7376/3 "green spec" blades*
- *Fully autoclaveable*
- *Guaranteed for life*

TimescoTM
SURGICAL + MEDICAL
INSTRUMENTS FOR LIFETM

THE NEW OPTIMA XL XENON HANDLE

See the Light

The Optima XL handle has been developed with one principle in mind - maximum illumination. The bright, white light is provided by a powerful new Xenon bulb. The unique new metal collar provides parabolic reflection of the light to ensure the maximum light emitted from the bulb is directed vertically upwards towards the fibre light guide in the blade. This ensures a minimum loss of light and so increases the light transmitted from the end of the fibre optic blade. Simple, but very effective.

New, easy bulb fitting

The bulb fitting has been radically redesigned to allow the bulb to be simply and easily changed, without the need for a bulb capsule. As the bulb now makes direct contact with the battery compartment, there is a minimal loss of current. This allows the majority of the battery power to be converted into light from the bulb. Unlike other systems,

there are no special tools required to replace the bulb.

Fully compatible

As the handle is made to the internationally recognised ISO 7376/3 standard, it may be used with not only our Optima fibre optic reusable blades, but also with our Freeway disposable plastic blades. With over 100% increase in light output, the combination of our Optima XL handle and our Freeway blades, is ideal for controlling cross infection, without compromising on the illumination of the airway and clinical performance.

Fully autoclaveable

Due to the increasing needs for decontamination to control cross infection, all our reusable laryngoscope blades and handles are fully autoclaveable. The new Optima XL handle is no exception.

'Lifetime Guarantee'

All Optima XL handles are produced from the highest quality materials and to the highest quality. Peace of mind comes as standard with our unrivalled Lifetime Guarantee. If the handle fails for any other reason than the bulb, we will repair or replace it free of charge, for life.



Why pay more?

Despite the fact that our Optima XL handles produce over 100% more light than similar handles, they are more competitively priced than other recognised laryngoscope brands. In some cases, almost half the price.

Try one yourself

Demonstrations and samples of these excellent new handles are available by contacting our Hospital Customer Support Team on 020 7511 1234 or from your local Territory Manager. Alternatively, click on our website www.timesco.com and we would be happy to arrange a sample for you.

Timesco Laryngoscopes - the best just keeps getting better.



3000.300.05
Optima XL Handle
Small (AA)

3000.300.10
Optima XL Handle
Medium (C)

3000.900.20
Optima XL Xenon Bulbs
Box of 6



Details on our complete range of Laryngoscope products may be found in our Laryngoscope Systems Catalogue. For samples or further information on any product, please contact our sales team on (020) 7511 1234.



Timesco
SURGICAL + MEDICAL
INSTRUMENTS FOR LIFE



Timesco of London, Timesco House, 1 Knights Road, London E16 2AT.
Telephone: +44 (020) 7511 1234 Fax: +44 (020) 7511 7888 Email: laryng@timesco.com

www.timesco.com

E & O.E. Copyright October 2002. Ref L.L.007/1
Optima and 'Instruments for Life' are registered trademarks of Timesco of London. All rights reserved.